



More information on the website
mirror.radwag.com/it/info,w1,JQP

Cilindrico, set (1 mg - 1 kg), imballaggi in plastica

OK-502-0018



The drawings, photos and graphics used are for illustrative purposes only.

Specifiche

Parametri metrologici	
Classe OIML	E1
Costruzione	
Materiale	Stainless steel
Peso	2 111,11 g
Valore nominale	1 mg – 1 kg
Numero	25 pcs

Compatibile con (Additional Fee)

Pinzette

Pennelli