



More information on the website  
[mirror.radwag.com/en/info,w1,JTA](https://mirror.radwag.com/en/info,w1,JTA)

# Stainless Steel Wire, set (1 mg - 500 mg), Wooden Box

OK-201-0041



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	F1
Construction	
Material	Stainless steel
Weight	1.11 g
Nominal value	1 mg - 500 mg
Weights quantity	12

## Compatible with (Additional Fee)

Tweezers

Brushes