



# Tips OMNITIP™ - Starter ST 1000 µl

PT-116-0001

More information on the website  
[mirror.radwag.com/en/info,w1,K9N](http://mirror.radwag.com/en/info,w1,K9N)



The drawings, photos and graphics used are for illustrative purposes only.

## Overview

OMNITIP™ Starter ST 1000 µl tips are high quality laboratory accessories designed for precision pipetting. Designed with efficiency and convenience in mind, they are an indispensable tool for everyday laboratory work.

Capacity: 1000 µl Colour: blue Set contains: box, 3 racks and 288 tips

Application: - precise dispensing of biological samples, - accurate pipetting of medicines and medical solutions, - preparing and transferring chemical solutions, - an indispensable tool in experiments requiring precise dosing.

## Compatible with (Additional Fee)

