



More information on the website
mirror.radwag.com/en/info,w1,KRZ

SAL/C Laboratory Antivibration Table



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Construction	
Construction	mild steel St3S
Load-bearing capacity	40 kg
Physical parameters	
Packaging dimensions	820×711×816 mm
Net weight	61 kg
Gross weight	72.5 kg
Length	770 mm
Width	680 mm
Height	815±10 mm

Compatible with

MA X7.IC.A Moisture Analyzer
 MUYA Modules
 MAS1 Modules
 MA 5Y Moisture Analyzer

MPS Modules
 C32.PM Precision Balance
 MA X7.A Moisture Analyzer
 MA 5Y.IC Moisture Analyzer

Device dimensions

SAL/C

