



More information on the website  
[mirror.radwag.com/en/info,w1,KRZ](http://mirror.radwag.com/en/info,w1,KRZ)

# SAL/C Laboratory Antivibration Table

WX-001-0076



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Construction	
Construction	mild steel St3S
Load-bearing capacity	40 kg
Physical parameters	
Packaging dimensions W x D x H	820×711×816 mm
Net weight	61 kg
Gross weight	72.5 kg
Length	770 mm
Width	680 mm
Height	815±10 mm

## Compatible with (Additional Fee)

MA X7.IC.A Moisture Analyzer  
 MA 5Y Moisture Analyzer  
 AS X7 Analytical Balance  
 PS R2 Precision Balance

MA R Moisture Analyzer  
 PS R2.H Precision Balance  
 XA 5Y.F Analytical Balance for Filters  
 PS X7 Precision Balance

WLC X7 Precision Balance  
C32.PM Precision Balance  
MA X7.A Moisture Analyzer  
MA 5Y.IC Moisture Analyzer

PS R1 Precision Balance  
AS R2 PLUS Analytical Balance  
AS R1 PLUS Analytical Balance  
MAS.1.Y Modules

## Device dimensions W x D x H

