



More information on the website  
[mirror.radwag.com/en/info,w1,KRZ](http://mirror.radwag.com/en/info,w1,KRZ)

# SAL/C Laboratory Antivibration Table

WX-001-0076



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Construction	
Construction	mild steel St3S
Load-bearing capacity	40 kg
Physical parameters	
Packaging dimensions	820×711×816 mm
Net weight	61 kg
Gross weight	72.5 kg
Length	770 mm
Width	680 mm
Height	815±10 mm

## Compatible with (Additional Fee)

MA X7.IC.A Moisture Analyzer  
 MUYA Modules  
 MAS1 Modules  
 MA 5Y Moisture Analyzer

MAS Modules  
 MPS Modules  
 C32.PM Precision Balance  
 MA X7.A Moisture Analyzer

MA X2.IC.A Moisture Analyzers  
AS X7 Analytical Balance  
WLC X7 Precision Balance  
MA X2.A Moisture Analyzers

PS X7 Precision Balance  
MA 5Y.IC Moisture Analyzer  
MA R Moisture Analyzer  
XA 5Y.F Analytical Balance for Filters

## Device dimensions

SAL/C

