



More information on the website
mirror.radwag.com/en/info,w1,KRZ

SAL/C Laboratory Antivibration Table

WX-001-0076



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Construction	
Construction	mild steel St3S
Load-bearing capacity	40 kg
Physical parameters	
Packaging dimensions	820×711×816 mm
Net weight	61 kg
Gross weight	72.5 kg
Length	770 mm
Width	680 mm
Height	815±10 mm

Compatible with (Additional Fee)

MA X7.IC.A Moisture Analyzer
 MUYA Modules
 MAS1 Modules
 MA 5Y Moisture Analyzer

MAS Modules
 MPS Modules
 C32.PM Precision Balance
 MA X7.A Moisture Analyzer

MA X2.IC.A Moisture Analyzers
AS X7 Analytical Balance
WLC X7 Precision Balance
MA X2.A Moisture Analyzers

PS X7 Precision Balance
MA 5Y.IC Moisture Analyzer
MA R Moisture Analyzer
XA 5Y.F Analytical Balance for Filters

Device dimensions

SAL/C

