



More information on the website
mirror.radwag.com/en/info,w1,KRZ

SAL/C Laboratory Antivibration Table

WX-001-0076



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Construction	
Construction	mild steel St3S
Load-bearing capacity	40 kg
Physical parameters	
Packaging dimensions W x D x H	820×711×816 mm
Net weight	61 kg
Gross weight	72.5 kg
Length	770 mm
Width	680 mm
Height	815±10 mm

Compatible with (Additional Fee)

MA X7.IC.A Moisture Analyzer
 MA 5Y Moisture Analyzer
 AS X7 Analytical Balance
 PS R2 Precision Balance

MA R Moisture Analyzer
 PS R2.H Precision Balance
 XA 5Y.F Analytical Balance for Filters
 PS X7 Precision Balance

WLC X7 Precision Balance
C32.PM Precision Balance
MA X7.A Moisture Analyzer
MA 5Y.IC Moisture Analyzer

PS R1 Precision Balance
AS R2 PLUS Analytical Balance
AS R1 PLUS Analytical Balance
MAS.1.Y Modules

Device dimensions W x D x H

