



More information on the website  
[mirror.radwag.com/en/info,w1,KV3](https://mirror.radwag.com/en/info,w1,KV3)

# Knob weights, Set (1 mg - 5 kg), Plastic Box

OK-102-0045



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	E2
Construction	
Material	Stainless steel
Weight	11 111.11 g
Nominal value	1 mg – 5 kg
Quantity	28 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes