



More information on the website
mirror.radwag.com/it/info,w1,KV3

Cilindrico, set (1 mg - 5 kg), imballaggi in plastica

OK-102-0045



The drawings, photos and graphics used are for illustrative purposes only.

Specifiche

Parametri metrologici	
Classe OIML	E2
Costruzione	
Materiale	Stainless steel
Peso	11 111,11 g
Valore nominale	1 mg – 5 kg
Numero	28 pcs

Compatibile con (Additional Fee)

Pinzette

Pennelli