



More information on the website
mirror.radwag.com/it/info,w1,LDM

Cilindrico, set (1 g - 2 kg), imballaggi in plastica

OK-102-0037



The drawings, photos and graphics used are for illustrative purposes only.

Specifiche

| Parametri metrologici | |
|-----------------------|-----------------|
| Classe OIML | E2 |
| Costruzione | |
| Materiale | Stainless steel |
| Peso | 6 110 g |
| Valore nominale | 1 g – 2 kg |
| Numero | 15 pcs |

Compatibile con (Additional Fee)

Pinzette

Pennelli