



More information on the website
mirror.radwag.com/de/info,w1,LDR

Prüfgewichte der M1 Klasse - Knopfgewichte, Set (1 g - 10 kg), Kunststoff-Etui

OK-402-0070



The drawings, photos and graphics used are for illustrative purposes only.

Technische Daten

Messtechnische Parameter	
OIML-Klasse	M1
Gewicht	21 110 g
Nominalwert	1 g – 10 kg
Anzahl	17 St.

Kompatibel mit (Additional Fee)

Pinsel