



More information on the website
mirror.radwag.com/en/info,w1,LDR

M1 Mass Standard, Knob Weights, Set (1 g - 10 kg), Plastic Box

OK-402-0070



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Metrological parameters

| | |
|---------------|-------------|
| OIML Class | M1 |
| Weight | 21 110 g |
| Nominal value | 1 g – 10 kg |
| Quantity | 17 pcs |

Compatible with (Additional Fee)

Brushes