



More information on the website  
[mirror.radwag.com/it/info,w1,LFY](https://mirror.radwag.com/it/info,w1,LFY)

# Peso standard F1 - cilindrico con cavità di adiacenza - 20 kg

OP-203-0097



The drawings, photos and graphics used are for illustrative purposes only.

## Specifiche

### Parametri metrologici

|             |    |
|-------------|----|
| Classe OIML | F1 |
|-------------|----|

### Costruzione

|           |                 |
|-----------|-----------------|
| Materiale | Stainless steel |
|-----------|-----------------|

|                 |       |
|-----------------|-------|
| Valore nominale | 20 kg |
|-----------------|-------|

## Compatibile con (Additional Fee)

Pennelli