



More information on the website  
[mirror.radwag.com/en/info,w1,LOG](https://mirror.radwag.com/en/info,w1,LOG)

# F2 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 mg - 500 g), Plastic Box



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

| Metrological parameters |                 |
|-------------------------|-----------------|
| OIML Class              | F2              |
| Construction            |                 |
| Material                | Stainless steel |
| Weight                  | 1 111.11 g      |
| Nominal value           | 1 mg - 500 g    |
| Quantity                | 24 pcs          |

## Compatible with

Tweezers

Brushes