

KIT 195.M Density Determination Kit

WX-002-0115





The drawings, photos and graphics used are for illustrative purposes only.

Overview

KIT 195.M is a density determination kit intended for PS M precisionbalances of high capacity, operating on the basis of RADWAG MonoBLOCK® technology, and featuring 195x195 weighing pan.


This stainless steel accessory allows determining density of:

solids, liquids.


Compatible with (Additional Fee)

PS X7 Balances for PGC PS R2 Precision Balance PS X7 Precision Balance PS 5Y Balances for PGC PS R1 Precision Balance PS 5Y Precision Balance


Device dimensions W x D x H









KIT 128



KIT 195



KIT 195.M



Granulte density determination KIT