



More information on the website  
[mirror.radwag.com/us/info,w1,LRO](https://mirror.radwag.com/us/info,w1,LRO)

# KIT 195.M Density Determination Kit

WX-002-0115



The drawings, photos and graphics used are for illustrative purposes only.

## Overview


Density determination kit is intended for PS M high capacity precision balances based on RADWAG MonoBLOCK technology equipped with 195×195 mm weighing pan. This stainless steel kit is applicable for determining density of: - solids; - liquids.

## Compatible with (Additional Fee)


PS X7 Balances for PGC  
 PS R2 Precision Balance  
 PS X7 Precision Balance

PS 5Y Balances for PGC  
 PS R1 Precision Balance  
 PS 5Y Precision Balance


Device dimensions W x D x H




KIT 85




KIT 128



KIT 195



KIT 195.M



Granulte density  
determination KIT