



More information on the website  
[mirror.radwag.com/en/info,w1,M3X](http://mirror.radwag.com/en/info,w1,M3X)

# F2 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 mg - 2 kg), Plastic Box

OK-302-0039



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	F2
Construction	
Material	Stainless steel
Weight	6 111.11 g
Nominal value	1 mg – 2 kg
Quantity	27 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes