



More information on the website  
[mirror.radwag.com/de/info,w1,MFX](https://mirror.radwag.com/de/info,w1,MFX)

# Prüfgewicht der M1 Klasse - Knopfgewicht - 5 kg

OP-403-0059



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

Messtechnische Parameter	
OIML-Klasse	M1
Konstruktion	
Material	Edelstahl
Nominalwert	5 kg

## Kompatibel mit (Additional Fee)

Pinsel