



More information on the website  
[mirror.radwag.com/tr/info,w1,MMB](http://mirror.radwag.com/tr/info,w1,MMB)

# Teraziler için Çanta PMV PLUS Microwave Moisture Analyzer

WX-003-0410



The drawings, photos and graphics used are for illustrative purposes only.

## Teknik Veriler

### Fiziksel parametreler

Ambalaj boyutları G x D x Y

760×510×425 mm

## İle uyumlu (Additional Fee)

PMV 5Y Nem Tayin Cihazı

## Cihaz boyutları G x D x Y

