



More information on the website  
[mirror.radwag.com/de/info,w1,MZS](http://mirror.radwag.com/de/info,w1,MZS)

# Prüfgewicht der M1 Klasse - Knopfgewicht - 200 g

OP-403-0125



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

| Messtechnische Parameter |           |
|--------------------------|-----------|
| OIML-Klasse              | M1        |
| Konstruktion             |           |
| Material                 | Edelstahl |
| Nominalwert              | 200 g     |

## Kompatibel mit (Additional Fee)

Pinzetten

Pinself