



Tips OMNITIP™ - Starter ST 300 µl

PT-115-0001

More information on the website
mirror.radwag.com/en/info,w1,N2A



The drawings, photos and graphics used are for illustrative purposes only.

Overview

OMNITIP™ Starter ST 300 µl tips are high quality laboratory accessories designed for precision pipetting. Designed with efficiency and convenience in mind, they are an indispensable tool for everyday laboratory work.

Capacity: 300 µl Colour: transparent Set contains: box, 6 racks and 576 tips

Application: - precise dispensing of biological samples, - accurate pipetting of medicines and medical solutions, - preparing and transferring chemical solutions, - an indispensable tool in experiments requiring precise dosing.

Compatible with (Additional Fee)

