



# Tips OMNITIP™ - Starter ST 300 µl

PT-115-0001

More information on the website  
[mirror.radwag.com/us/info,w1,N2A](https://mirror.radwag.com/us/info,w1,N2A)



The drawings, photos and graphics used are for illustrative purposes only.

## Overview

OMNITIP™ - Starter ST 300,transparent, 300 µl,box + 6 racks + 576 tips.

## Compatible with (Additional Fee)

Automatic Variable-Volume Pipettes

Workstation for pipettes calibration