



More information on the website
mirror.radwag.com/de/info,w1,N88

Prüfgewichte der M1 Klasse - Knopfgewichte, Set (1 mg - 10 kg), Holz- Etui

OK-401-0086



The drawings, photos and graphics used are for illustrative purposes only.

Technische Daten

Messtechnische Parameter	
OIML-Klasse	M1
Gewicht	21 111,11 g
Nominalwert	1 mg – 10 kg
Anzahl	29 St.

Kompatibel mit (Additional Fee)

Pinsel