



More information on the website  
[mirror.radwag.com/it/info,w1,NBC](https://mirror.radwag.com/it/info,w1,NBC)

# Peso standard F1 - filo di acciaio inossidabile - 20 mg

OP-203-0070



The drawings, photos and graphics used are for illustrative purposes only.

## Specifiche

Parametri metrologici	
Classe OIML	F1
Costruzione	
Materiale	Stainless steel
Valore nominale	20 mg

## Compatibile con (Additional Fee)

Pinzette

Pennelli