



More information on the website  
[mirror.radwag.com/en/info,w1,NFX](https://mirror.radwag.com/en/info,w1,NFX)

# F1 Mass Standard - Knob Weights Without Adjustment Chamber, Set (1 g - 500 g), Wooden Box

OK-201-0064



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	F1
Construction	
Material	Stainless steel
Weight	1 110 g
Nominal value	1 – 500 g
Quantity	12 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes