



More information on the website
mirror.radwag.com/en/info,w1,NNM

Under-pan Weighing Rack for XA Balances

WX-002-0003



The drawings, photos and graphics used are for illustrative purposes only.

Overview

A rack for under-pan weighing intended for XA balances.

The rack is designed for weighing loads requiring carrying out the measuring process outside the main weighing platform, e.g. loads with non-standard dimensions and sizes and those generating a magnetic field.

The rack also facilitates determining density of loads.

The rack is available in powder coated mild steel finishing.

Compatible with (Additional Fee)

XA 5Y.A Analytical Balance

XA 5Y.F Analytical Balance for Filters

XA 5Y Analytical Balance