



# Under-pan Weighing Rack for XA Balances

WX-002-0003

More information on the website  
[mirror.radwag.com/en/info,w1,NNM](https://mirror.radwag.com/en/info,w1,NNM)



The drawings, photos and graphics used are for illustrative purposes only.

## Overview

A rack for under-pan weighing intended for XA balances.

The rack is designed for weighing loads requiring carrying out the measuring process outside the main weighing platform, e.g. loads with non-standard dimensions and sizes and those generating a magnetic field.

The rack also facilitates determining density of loads.

The rack is available in powder coated mild steel finishing.

## Compatible with (Additional Fee)

XA 5Y.A Analytical Balance  
XA 5Y.F Analytical Balance for Filters

XA 5Y Analytical Balance