



More information on the website
mirror.radwag.com/us/info,w1,NNM

Under-hook weighing rack for XA balances

WX-002-0003



The drawings, photos and graphics used are for illustrative purposes only.

Overview

A rack for under-hook weighing intended for balances XA.

The rack is designed for weighing loads requiring carrying out the measuring process outside the main weighing platform, e.g. loads with non-standard dimensions and sizes and those generating a magnetic field.

The rack also facilitates determining density of loads.

The rack is available in powder coated mild steel finishing.

Compatible with (Additional Fee)

XA 5Y.A Analytical Balance
Analytical Balance for Filters XA 5Y.F

XA 5Y Analytical Balance