



More information on the website  
[mirror.radwag.com/de/info,w1,NUT](https://mirror.radwag.com/de/info,w1,NUT)

# Windschutzschrank für XA 4Y und XA 5Y Waagen

WX-001-0160



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

### Physikalische Parameter


Dimension	560x350x352 mm
Verpackungsgröße B x T x H	695x485x485 mm
Nettogewicht	11 kg
Bruttogewicht	16 kg

## Kompatibel mit (Additional Fee)

Mikrowaage für Pipettenkalibrierung XA 5Y.M.A.P  
 Mikrowaage XA 5Y.M.A  
 Analysenwaage XA 5Y.A  
 Analysenwaage XA 5Y


Manueller Massekomparator XA 5Y.KO  
 Mikrowaage XA 5Y.M  
 Mikrowaage für Stentwägung XA 5Y.M.A.S

# Abmessungen des Geräts B x T x H




Anti-draft chamber for balances with 128x128 mm weighing pan


Anti-draft chamber for X2 and R2 precision balances




Anti-draft chamber for microbalances




Draft shield for precision balances




Anti-draft chamber for APP-00-SK mass comparators



Anti-draft chamber for APP-KO-SK mass comparators



Anti-draft chamber for WAY-KO mass comparators



Anti-draft chamber for XA 5Y / XA 5Y.A balances