



More information on the website
mirror.radwag.com/en/info,w1,NUT

Anti-Draft Chamber for XA 4Y and XA 5Y Balances

WX-001-0160



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Physical parameters

Dimension	560x350x352 mm
Packaging dimensions W x D x H	695x485x485 mm
Net weight	11 kg
Gross weight	16 kg

Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration
 XA 5Y.M.A Microbalance
 XA 5Y.A Analytical Balance
 XA 5Y Analytical Balance

XA 5Y.KO Manual Mass Comparator
 XA 5Y.M Microbalance
 XA 5Y.M.A.S Microbalance for stents weighing

Device dimensions W x D x H



Anti-draft chamber for balances with 128x128 mm weighing pan

Anti-draft chamber for X2 and R2 precision balances



Anti-draft chamber for microbalances



Draft shield for precision balances



Anti-draft chamber for APP-00-SK mass comparators



Anti-draft chamber for APP-KO-SK mass comparators



Anti-draft chamber for WAY-KO mass comparators



Anti-draft chamber for XA 5Y / XA 5Y.A balances