



More information on the website  
[mirror.radwag.com/en/info,w1,NUT](https://mirror.radwag.com/en/info,w1,NUT)

# Anti-Draft Chamber for XA 4Y and XA 5Y Balances



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Physical parameters	
Dimension	560x350x352 mm
Packaging dimensions	695x485x485 mm
Net weight	11 kg
Gross weight	16 kg

## Compatible with

XA 5Y.M.A.P Microbalance for pipettes calibration  
 XA 5Y.M.A Microbalance  
 XA 5Y.A Analytical Balance  
 XA 5Y Analytical Balance

XA 5Y.KO Manual Mass Comparator  
 XA 5Y.M Microbalance  
 XA 5Y.M.A.S Microbalance for stents weighing

## Device dimensions




Anti-draft chamber for balances with 128x128 mm weighing pan


Anti-draft chamber for X2 and R2 precision balances



Anti-draft chamber for microbalances



Draft shield for precision balances



Anti-draft chamber for APP-00-SK mass comparators



Anti-draft chamber for APP-KO-SK mass comparators



Anti-draft chamber for WAY-KO mass comparators



Anti-draft chamber for XA 5Y / XA 5Y.A balances