



More information on the website  
[mirror.radwag.com/fr/info,w1,NUT](https://mirror.radwag.com/fr/info,w1,NUT)

# Boîte anti-courant d'air pour XA 4Y/ XA 5Y balances

WX-001-0160



The drawings, photos and graphics used are for illustrative purposes only.

## Paramètres Techniques

### Paramètres physiques

|                               |                |
|-------------------------------|----------------|
| Dimension                     | 560x350x352 mm |
| Dimensions de colis L x P x H | 695x485x485 mm |
| Masse nette                   | 11 kg          |
| Masse brute                   | 16 kg          |

## Compatible avec (Additional Fee)

Microbalance pour étalonnage des pipettes XA 5Y.M.A.P  
 Microbalance XA 5Y.M.A  
 Balance analytique XA 5Y.A  
 Balance analytique XA 5Y

Comparteur manuel de masse XXA 5Y.KO  
 Microbalance XA 5Y.M  
 Microbalance pour stents XA 5Y.M.A.S

# Dimensions d'appareil L x P x H



Anti-draft chamber for balances with 128x128 mm weighing pan

Anti-draft chamber for X2 and R2 precision balances



Anti-draft chamber for microbalances



Draft shield for precision balances



Anti-draft chamber for APP-00-SK mass comparators



Anti-draft chamber for APP-KO-SK mass comparators



Anti-draft chamber for WAY-KO mass comparators



Anti-draft chamber for XA 5Y / XA 5Y.A balances