



More information on the website
mirror.radwag.com/fr/info,w1,NUT

Boîte anti-courant d'air pour XA 4Y/ XA 5Y balances

WX-001-0160



The drawings, photos and graphics used are for illustrative purposes only.

Paramètres Techniques

Paramètres physiques

Dimension	560x350x352 mm
Dimensions de colis	695x485x485 mm
Masse nette	11 kg
Masse brute	16 kg

Compatible avec (Additional Fee)

Microbalance pour étalonnage des pipettes XA 5Y.M.A.P
 Microbalance XA 5Y.M.A
 Balance analytique XA 5Y.A
 Balance analytique XA 5Y

Comparteur manuel de masse XXA 5Y.KO
 Microbalance XA 5Y.M
 Microbalance pour stents XA 5Y.M.A.S

Dimensions d'appareil



Anti-draft chamber for balances with 128x128 mm weighing pan

Anti-draft chamber for X2 and R2 precision balances




Anti-draft chamber for microbalances




Draft shield for precision balances




Anti-draft chamber for APP-00-SK mass comparators



Anti-draft chamber for APP-KO-SK mass comparators



Anti-draft chamber for WAY-KO mass comparators



Anti-draft chamber for XA 5Y / XA 5Y.A balances