



More information on the website  
[mirror.radwag.com/it/info,w1,NUT](https://mirror.radwag.com/it/info,w1,NUT)

# Anti-Draft Chamber for XA 4Y and XA 5Y Balances

WX-001-0160



The drawings, photos and graphics used are for illustrative purposes only.

## Specifiche

### Parametri fisici

Dimensioni	560x350x352 mm
Dimensioni del pacco	695x485x485 mm
Peso net	11 kg
Peso lordo	16 kg

## Compatibile con (Additional Fee)

Microbilancia per la calibrazione delle pipette XA 5Y.M.A.P  
 Microbilancia XA 5Y.M.A  
 Bilancia analitica XA 5Y.A  
 Bilancia analitica XA 5Y

XA 5Y.KO Comparatore di massa manuale  
 Microbilancia XA 5Y.M  
 Microbilancia per stent XA 5Y.M.A.S

## Device dimensions




Anti-draft chamber for balances with 128x128 mm weighing pan


Anti-draft chamber for X2 and R2 precision balances




Anti-draft chamber for microbalances




Draft shield for precision balances




Anti-draft chamber for APP-00-SK mass comparators



Anti-draft chamber for APP-KO-SK mass comparators



Anti-draft chamber for WAY-KO mass comparators



Anti-draft chamber for XA 5Y / XA 5Y.A balances