



More information on the website
mirror.radwag.com/it/info,w1,NUT

Anti-Draft Chamber for XA 4Y and XA 5Y Balances

WX-001-0160



The drawings, photos and graphics used are for illustrative purposes only.

Specifiche

Parametri fisici

Dimensioni	560x350x352 mm
Dimensioni del pacco L x P x A	695x485x485 mm
Peso net	11 kg
Peso lordo	16 kg

Compatibile con (Additional Fee)

Microbilancia per la calibrazione delle pipette XA 5Y.M.A.P
 Microbilancia XA 5Y.M.A
 Bilancia analitica XA 5Y.A
 Bilancia analitica XA 5Y

XA 5Y.KO Comparatore di massa manuale
 Microbilancia XA 5Y.M
 Microbilancia per stent XA 5Y.M.A.S

Device dimensions L x P x A



Anti-draft chamber for balances with 128x128 mm weighing pan

Anti-draft chamber for X2 and R2 precision balances



Anti-draft chamber for microbalances



Draft shield for precision balances




Anti-draft chamber for APP-00-SK mass comparators



Anti-draft chamber for APP-KO-SK mass comparators



Anti-draft chamber for WAY-KO mass comparators



Anti-draft chamber for XA 5Y / XA 5Y.A balances