



More information on the website
mirror.radwag.com/tr/info,w1,NUT

Anti-Draft Chamber for XA 4Y and XA 5Y Balances

WX-001-0160



The drawings, photos and graphics used are for illustrative purposes only.

Teknik Veriler

Fiziksel parametreler


Boyut	560x350x352 mm
Ambalaj boyutları G x D x Y	695x485x485 mm
Net ağırlık	11 kg
Brüt ağırlık	16 kg

İle uyumlu (Additional Fee)

XA 5Y.M.A.P pipetlerini kalibre etmek için Mikro Terazi
XA 5Y.M.A Mikro Teraziler
XA 5Y.A Analitik Terazi
XA 5Y Analitik Terazi


XA 5Y.KO Manuel kütle komparatör
XA 5Y.M Mikro Teraziler
XA 5Y.M.A.S Stentler için Mikro Teraziler

Cihaz boyutları G x D x Y



Anti-draft chamber for balances with 128x128 mm weighing pan

Anti-draft chamber for X2 and R2 precision balances




Anti-draft chamber for microbalances




Draft shield for precision balances




Anti-draft chamber for APP-00-SK mass comparators



Anti-draft chamber for APP-KO-SK mass comparators



Anti-draft chamber for WAY-KO mass comparators



Anti-draft chamber for XA 5Y / XA 5Y.A balances