



More information on the website
mirror.radwag.com/us/info,w1,NUT

Anti-Draft Chamber for XA 4Y and XA 5Y Balances

WX-001-0160



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Physical parameters


| | |
|--------------------------------|----------------|
| Dimension | 560x350x352 mm |
| Packaging dimensions W x D x H | 695x485x485 mm |
| Net weight | 11 kg |
| Gross weight | 16 kg |

Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration
 XA 5Y.M.A Microbalance
 XA 5Y.A Analytical Balance
 XA 5Y Analytical Balance


XA 5Y.KO Manual Mass Comparator
 XA 5Y.M Microbalance
 XA 5Y.M.A.S Microbalance for stents weighing

Device dimensions W x D x H




Anti-draft chamber for balances with 128x128 mm weighing pan


Anti-draft chamber for X2 and R2 precision balances




Anti-draft chamber for microbalances




Draft shield for precision balances




Anti-draft chamber for APP-00-SK mass comparators



Anti-draft chamber for APP-KO-SK mass comparators



Anti-draft chamber for WAY-KO mass comparators



Anti-draft chamber for XA 5Y / XA 5Y.A balances