



More information on the website
mirror.radwag.com/us/info,w1,ODT

DJ-04 Antistatic ionizer

WX-001-0097



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Power supply	100 – 240 V AC 50/60 Hz
Operating temperature	0 – +40 °C
Relative humidity	5% – 80%
Distance form a sample	Ok. 5 – 50 mm
Ion source	stainless steel
Physical parameters	
Frame dimensions	268×110×295 mm
Packaging dimensions W x D x H	375×375×130 mm
Net weight	1,3 kg
Gross weight	2,5 kg

Compatible with (Additional Fee)


XA 5Y.M.A.P Microbalance for pipettes calibration
 AS X7 Analytical Balance

XA 5Y Analytical Balance
 AS R2 PLUS Analytical Balance


XA 5Y.M.A Microbalance
XA 5Y.A Analytical Balance
AS 5Y Analytical Balance
Analytical Balance for Filters XA 5Y.F

XA 5Y.M Microbalance
AS R1 PLUS Analytical Balance
XA 5Y.M.A.S Microbalance for stents weighing


Device dimensions W x D x H



DJ-04



DJ-05



power supply