



More information on the website
mirror.radwag.com/de/info,w1,OGP

Pumpe zum Absaugen von Flüssigkeiten für AP-12/12.5Y

WX-001-0480



The drawings, photos and graphics used are for illustrative purposes only.

Beschreibung

Dedicated Pump Set The pump enables simultaneous liquid suction from all measuring channels. As a result, all vessels can be emptied at the same time, without the need to perform this operation separately for each one. This solution significantly improves workflow and shortens measurement time when calibrating multichannel pipettes. Simultaneous liquid suction also increases user convenience and improves the overall flow of the calibration process.

The suction pump for liquids is a standard (free) component of the AP-12/12.5Y Automatic Device for Multichannel Pipette Calibration.

Kompatibel mit (Additional Fee)

Automatische Waagen für die Pipettenkalibrierung

