



## Suction pump for liquids for AP-12/12.5Y

WX-001-0480

More information on the website  
[mirror.radwag.com/en/info,w1,OGP](http://mirror.radwag.com/en/info,w1,OGP)



The drawings, photos and graphics used are for illustrative purposes only.

### Overview

**Dedicated Pump Set** The pump enables simultaneous liquid suction from all measuring channels. As a result, all vessels can be emptied at the same time, without the need to perform this operation separately for each one. This solution significantly improves workflow and shortens measurement time when calibrating multichannel pipettes. Simultaneous liquid suction also increases user convenience and improves the overall flow of the calibration process.

The suction pump for liquids is a standard (free) component of the AP-12/12.5Y Automatic Device for Multichannel Pipette Calibration.

### Compatible with (Additional Fee)

Automatic Balances for Pipettes Calibration