



More information on the website  
[mirror.radwag.com/en/info,w1,0SJ](https://mirror.radwag.com/en/info,w1,0SJ)

# F1 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 kg - 5 kg), Plastic Box

OK-202-0062



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	F1
Construction	
Material	Stainless steel
Weight	10 000 g
Nominal value	1 – 5 kg
Quantity	4 pcs

## Compatible with (Additional Fee)

Brushes