



More information on the website  
[mirror.radwag.com/en/info,w1,PDJ](http://mirror.radwag.com/en/info,w1,PDJ)

# M1 Mass Standard, Knob Weights, Set (1 mg - 500 g), Plastic Box

OK-402-0074



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

### Metrological parameters

OIML Class	M1
Weight	1 111.11 g
Nominal value	1 mg - 500 g
Quantity	24 pcs

## Compatible with (Additional Fee)

Brushes