



More information on the website  
[mirror.radwag.com/en/info,w1,PEX](http://mirror.radwag.com/en/info,w1,PEX)

# E1 Mass Standard - knob weight - 2 g

OP-503-0047



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

| Metrological parameters |                 |
|-------------------------|-----------------|
| OIML Class              | E1              |
| Construction            |                 |
| Material                | Stainless steel |
| Nominal value           | 2 g             |

## Compatible with (Additional Fee)

Tweezers

Brushes