



More information on the website
mirror.radwag.com/en/info,w1,PJV

F1 Mass Standard - Knob Weights Without Adjustment Chamber, Set (1 g - 2 kg), Wooden Box

OK-201-0066



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Metrological parameters	
OIML Class	F1
Construction	
Material	Stainless steel
Weight	6 110 g
Nominal value	1 g – 2 kg
Quantity	15 pcs

Compatible with (Additional Fee)

Tweezers

Brushes