



More information on the website  
[mirror.radwag.com/de/info,w1,PZ0](https://mirror.radwag.com/de/info,w1,PZ0)

# Prüfgewicht der M1 Klasse - Knopfgewicht - 10 g

OP-403-0121



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

Messtechnische Parameter	
OIML-Klasse	M1
Konstruktion	
Material	Edelstahl
Nominalwert	10 g

## Kompatibel mit (Additional Fee)

Pinzetten

Pinself