



# M1 Mass Standard - Knob Weights - 10 g

OP-403-0121

More information on the website  
[mirror.radwag.com/en/info,w1,PZ0](http://mirror.radwag.com/en/info,w1,PZ0)



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

| Metrological parameters |                 |
|-------------------------|-----------------|
| OIML Class              | M1              |
| Construction            |                 |
| Material                | Stainless steel |
| Nominal value           | 10 g            |

## Compatible with (Additional Fee)

Tweezers

Brushes