



More information on the website
mirror.radwag.com/en/info,w1,QKR

Stand for terminal 5Y / PUE CY10 - mild steel version

WX-004-0298



The drawings, photos and graphics used are for illustrative purposes only.

Overview

A stand for supporting scale's indicator installed on a cable 5Y / PUE CY10.
Mild steel version available.

Compatible with (Additional Fee)

5Y.PM Precision Balance
CY10 Scales for PGC
CY10 Multifunctional Scale

CY10 Precision Balance
HRP 5Y.KO Mass Comparator
HRP 5Y.KB Mass Comparator

Device dimensions W x D x H

