



More information on the website  
[mirror.radwag.com/us/info,w1,QKR](https://mirror.radwag.com/us/info,w1,QKR)

# Stand for terminal 5Y / PUE CY10 - mild steel version

WX-004-0298



The drawings, photos and graphics used are for illustrative purposes only.

## Overview

A stand, based on a tripod for supporting scale's indicator installed on a cable.

## Compatible with (Additional Fee)

5Y.PM Precision balance  
 CY10 Scales for PGC  
 CY10 Multifunctional Scale

CY10 Precision balance  
 HRP 5Y.KO Mass Comparator  
 HRP 5Y.KB Mass Comparator

## Device dimensions W x D x H

