

## Stand for terminal 5Y / PUE CY10 - mild steel version

WX-004-0298





The drawings, photos and graphics used are for illustrative purposes only.

## **Overview**

A stand, based on a tripod for supporting scale's indicator installed on a cable.

## **Compatible with (Additional Fee)**

5Y.PM Precision balance CY10 Scales for PGC CY10 Multifunctional Scale CY10 Precision balance HRP 5Y.KO Mass Comparator HRP 5Y.KB Mass Comparator

## Device dimensions W x D x H

